

TechLITE Pen Needles Compatibility Chart

For pen needles made in China.

TechLITE Pen Needles fit most standard insulin pens.

Standard Needle-based Injection System EN ISO 11608-2:2012 clause 4.9 and 11	
Sanofi-Aventis	ClickSTAR [®]
	SoloSTAR [®]
	JuniorSTAR [®]
Eli Lilly	KwikPen [®]
	HumaPen [®] Memoir
	HumaPen [®] LUXURA HD
	HumaPen [®] LUXURA
	HumaPen [®] SAVVIO
	HumaPen [®] Ergo II
Novo Nordisk	FlexTOUCH [®]
	NovoPen [®] 3
	NovoPen [®] 4
	NovoPen [®] 5
	Innolet [®]
	FlexPen [®]
	Victoza [®]
	NovoPen Echo [®]
Owen Mumford	Autopen [®] Classic, Autopen [®] 24
Berlin Chemie	BerliPen [®] Areo
	BerliPen [®] Junior

NOTE! The list above is not exhaustive - the brand names may vary in different markets, change due to mergers or acquisitions by their owners and private label arrangements. All specifications are subject to change. Consult appropriate company for additional information. Trademarks are the property of their respective owners.